

Infu-Disk™ Order Form

Qty	Vol(mL)	Rate (mL/hr)	Connector*	Sterile	Model #	Price (ea)	Total
			<input type="checkbox"/> Male Luer <input type="checkbox"/> Cannula	<input type="checkbox"/> Yes <input type="checkbox"/> No		\$41.00 All sterile, standard rates	
			<input type="checkbox"/> Male Luer <input type="checkbox"/> Cannula	<input type="checkbox"/> Yes <input type="checkbox"/> No		\$44.00 All sterile, non-standard rates and/or connectors	
			<input type="checkbox"/> Male Luer <input type="checkbox"/> Cannula	<input type="checkbox"/> Yes <input type="checkbox"/> No		\$39.00 All non-sterile, standard rates	
			<input type="checkbox"/> Male Luer <input type="checkbox"/> Cannula	<input type="checkbox"/> Yes <input type="checkbox"/> No		\$42.00 All non-sterile, non-standard rates and/or connectors	

*Male Luer is standard for all configurations

CA tax 7.5% (if applicable) _____

Total

(does not include shipping charges)

PO Number _____

Company _____

Contact _____

Phone () _____ - _____

Fax () _____ - _____

E-Mail _____

Shipping Address

Street _____

City _____

State/Prov _____

Zip/Postal Code _____

Country _____

Billing Address same

Street _____

City _____

State/Prov _____

Zip/Postal Code _____

Country _____

Shipping Method Fed Ex 2 Day Fed Ex Std Overnight USPS Express (1 day)
 USPS Priority (2 day) Fed Ex Express Saver (3 day)

Expected Arrival Date / /

4 to 6 weeks delivery for pumps not in stock. Expedite Fee (delivery in 10 working days) \$410 plus 10% per pump for orders up to 200 pumps. Contact Med-e-Cell for larger orders.

Comments _____

How to Order

Standard* Infu-Disk™ Configurations:

Model	Description
05-0030-S	5 mL, 0.03 mL/hr (7 day delivery)
05-0042-S	5 mL, 0.042 mL/hr (5 day delivery)
05-0070-S	5 mL, 0.07 mL/hr (3 day delivery)
05-0110-S	5 mL, 0.11 mL/hr (2 day delivery)
05-0210-S	5 mL, 0.21 mL/hr (1 day delivery)
10-0060-S	10 mL, 0.06 mL/hr (7 day delivery)
10-0083-S	10 mL, 0.083 mL/hr (5 day delivery)
10-0140-S	10 mL, 0.14 mL/hr (3 day delivery)
10-0210-S	10 mL, 0.21 mL/hr (2 day delivery)
10-0420-S	10 mL, 0.42 mL/hr (1 day delivery)

* flow rates are for 25°C with 6" H₂O fluid pressure

Custom Infu-Disk™ Configurations:



- ① 05 or 10 (volume in mL)
- ② Flow Rate in mL/hr x 1000
- ③ S (sterile) or N (non - sterile)

Example: Model 10-0350-S refers to a sterile 10 mL pump with a flow rate of 0.35 mL/hr.

5 mL pump flow rate range: 0.03 to 0.5 mL/hr
10 mL pump flow rate range: 0.03 to 4.0 mL/hr

Order Infu-Disk™ pumps by phone: **858-552-0781** fax: **858-552-8018** or online: **www.medecell.com**



10633 Roselle St., San Diego, CA 92121

Infu-Disk™ U.S. Patent #4,902,278 et al.